

# MobilityatSea

## Moving & Handling

Reference No. MH021

Cricket II Transfer Aid



This Transfer Aid is not suitable for transferring to an electric wheelchair unless there is clear access underneath (ground clearance of 12cm or above) – see specifications below.

The user must be able to stand up and remain in a standing/sitting position with adequate arm strength to use the crossbar.

### Specifications

Overall Length	89cm / 35"
Overall Height	100cm / 39.3"
Height of Base Plate	10cm / 3.9"
Height of Knee Pad (adjustable)	46 - 65cm / 18.1 – 25.5"
Width Between Legs (front)	24.5cm / 9.6"
Total Width Legs (front)	32cm / 12.6"
Width between arms (max)	52.5cm / 20.6"
Height of back support	92cm / 36.2"
Overall Width	68cm / 26.7"
Weight	25kg / 55.1lbs / 3.9st
Max User Weight	133.3kg / 294lbs / 21st

### **Disclaimer**

Certain circumstances may determine that another product is supplied. We will endeavour to notify you in such an event however this will not impact on the functionality of the product or the purpose of which is it intended.