



## Moving & Handling

Reference No. MH023

Drive Transfer Aid



This Transfer Aid has manual adjustable legs that splay open for closer transfers. **The user must be able to stand up and remain standing until the paddles are lowered into position and vice versa.**

### Specifications

<b>Overall Height</b>	105cm / 41.3"	<b>Height of Handle (from base)</b>	91cm / 35.8"
<b>Overall Length</b>	80cm / 31.4"	<b>Height of Seat (min)</b>	52cm / 20.4"
<b>Overall Width (min)</b>	52cm / 20.4"	<b>Height of Seat (max)</b>	71cm / 27.9"
<b>Overall Width (max)</b>	83cm / 32.6"	<b>Height of Knee Pad (adjustable)</b>	24 – 47cm / 9.4 – 18.5"
<b>Height of Base Plate</b>	10cm / 4"	<b>Width Between Front Legs (min)</b>	43cm / 17"
<b>Weight</b>	30.6kg / 67.4lbs / 4st 8	<b>Width Between Front Legs (max)</b>	76.2cm / 30"
<b>Max User Weight</b>	177.8kg / 392lbs / 28st	<b>Seat Width</b>	46cm / 18.1"

### **Disclaimer**

Certain circumstances may determine that another product is supplied. We will endeavour to notify you in such an event however this will not impact on the functionality of the product or the purpose of which is it intended.