

## Wheelchairs

Reference no. WC004

Heavy Duty Transit Wheelchair

(22" with brakes on the handles)



This wheelchair is designed for the larger passenger. The brakes are positioned on the handles for the attendant to control.

We advise that you plan your route to avoid adverse terrain, inclines, kerbs and steps (if this is unavoidable, we do not advise an incline above 5 degrees and kerb or step above 3cm/1.1").

**Please be mindful that gangways and airbridges can vary in the degree of gradient due to the tides, some products are not designed to cope with acute slopes.**

**THIS PRODUCT NEEDS TO BE FOLDED AND PUSHED THROUGH NARROW CABIN DOORWAYS.**

### Specifications

<b>Height</b>	93.9cm / 37"
<b>Width</b>	71.1cm / 28"
<b>Closed Width</b>	26.6cm / 10.5"
<b>Length</b>	102.5cm / 40.3"
<b>Seat Height</b>	48.2cm / 19"
<b>Seat Width</b>	55.8cm / 22"
<b>Seat Depth</b>	48.2cm / 19"
<b>Weight</b>	23kg / 51.5lbs / 3st 6
<b>Max User Weight</b>	196.8kg / 434lbs / 31st

### Disclaimer

Certain circumstances may determine that another product is supplied. We will endeavour to notify you in such an event however this will not impact on the functionality of the product or the purpose of which is it intended.