

Moving & Handling

Reference no. MH027

Molift Raiser Pro



The Molift Raiser Pro allows for a sit-to-stand transfer – **the user must be able to stand up and weight bear.**

Due to the thickness of the carpeting on board, it is not recommended to travel across the floor but to simply use it to turn the user 90 – 180°.

Specifications

Overall Height	124cm / 48.8"	Height Of Handle (top – from baseplate)	123cm / 48.4"
Overall Width	53cm / 20.8"	Height Of Handle (middle – from baseplate)	108cm / 42.5"
Overall Length	62cm / 24.4"	Height Of Handle (bottom – from baseplate)	92cm / 36.2"
Weight	10kg / 22lbs / 1.5st	Height Of Kneepad (adjustable)	32 - 53cm / 12.5 – 20.8"
Max User Weight	158.7kg / 350lbs / 25st	Width Between Arms	40.5cm / 15.9"

Disclaimer

Certain circumstances may determine that another product is supplied. We will endeavour to notify you in such an event however this will not impact on the functionality of the product or the purpose of which is it intended.